CBC PlasmaEgg®

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CE

REF 500610
As of 2024-04
GB500610GB

The logo for Kimetec, featuring the brand name in a bold, sans-serif font with a stylized leaf icon above the letter 'i'.

 **Please read the instructions for use carefully before using the PlasmaEgg® and store them with the device.**

Content

1.	Introduction and general information	19
1.1	Disinfection principle	19
2.	Intended use	19
2.1	Performance characteristics	19
2.2	Material compatibility	19
3.	General safety instructions	20
3.1	Operating company's duty of care	20
4.	Structure and components	21
4.1	Control panel	21
4.2	Device	21
5.	Transport and installation	22
5.1	Transport	22
5.2	Installation	22
6.	Disinfection with the PlasmaEgg®	22
6.1	Preparation	22
6.2	Switching on the device	22
6.3	Selecting a programme	22
6.4	Starting the programme	22
6.5	Removing the disinfected item	23
6.6	Switching off the device	23
7.	Replacing consumables	23
8.	Hygiene	24
8.1	Cleaning	24
8.2	Disinfection	24
9.	Conditions for storage and transport and ambient conditions	24
10.	Warranty	24
11.	Liability	24
12.	Troubleshooting (faults and remedies)	25
13.	Environmental protection and disposal	25
14.	Technical specifications	25
15.	Electromagnetic compatibility (EMC)	26
15.1	Behaviour of the PlasmaEgg® in the described electromagnetic environment	26
15.2	Description of the electromagnetic environment	26
15.3	Emissions	27
15.4	Immunity	27
16.	Legend	29
17.	Accessories and spare parts	30

1. Introduction and general information

1.1 Disinfection principle

The PlasmaEgg® disinfection device generates cold atmospheric plasma. Cold plasma generates ozone, which contributes to the disinfectant effect.

2. Intended use

- The PlasmaEgg® is intended exclusively for its intended use: the disinfection of everyday objects, laboratory equipment and personal protective equipment. Any other use is prohibited.
- The item being disinfected must be suitable for disinfection with cold plasma.

2.1 Performance characteristics

The following table shows the action spectra that are reached after the relevant treatment times:

Performance spectrum	
Testing and classification based on the specifications of DIN EN 14885, DIN EN 17111 and DIN EN 16777 for chemical disinfectants and antiseptics	
Disinfection of non-porous surfaces*	
Bactericidal	20 minutes
Yeasticidal	20 minutes
Fungicidal	60 minutes
Virucidal	90 minutes
Disinfection of porous surfaces**	
Bactericidal	20 minutes
Yeasticidal	20 minutes
Mycobactericidal	20 minutes
Fungicidal	90 minutes
Virucidal	20 minutes

* e.g. silicone, plastic (tested without interfering substances)

** e.g. fabric (tested without interfering substances)

The programme running time consists of the disinfection time (performance spectrum) plus a 5-minute plasma break-down phase.

2.2 Material compatibility

- Glass or corrosion-resistant metals are resistant to cold plasma.
- Untreated metals and certain plastics may experience surface corrosion or colour changes as a result of being treated with ozone.
- Tables of the ozone compatibility of materials can provide an indication of compatibility. We recommend testing to verify compatibility.

3. General safety instructions

- The device must not be opened during operation.
- If a fault occurs during disinfection, only open the lid in a ventilated room with a volume greater than 36 m³. Keep a minimum distance of 60 cm between your face and the device when opening it.
- The disinfection process was only successfully completed when the display shows “000”.
- The maximum ozone concentration generated in the chamber during disinfection is approximately 1200 ppm. Under normal operating conditions, the ozone concentration liberated from the chamber is below the permissible concentration at the workplace for 8 hour exposure.

If the operator opens the lid during disinfection, the average ozone concentration liberated from the chamber is 2.7 ppm. This corresponds to 0.1 ppm in a standard room (36 m³) which complies to the permissible ozone concentration at the workplace for short term exposure (0.1 ppm).

- The device must not be operated outside or while being transported.
- The device must not be operated in potentially explosive environments or oxygen-enriched atmospheres.
- The device must not be operated in the vicinity of flammable anaesthetics or other flammable substances.
- The device must not be operated if the housing or mains cable are damaged.
- The device must not be opened during operation.
- Do not spray or pour any liquid directly onto or into the device.
- Always disconnect the mains plug before cleaning and/or disinfecting the device.
- Only use accessories and spare parts that meet the manufacturer’s specifications.
- Keep the device away from children.

3.1 Operating company’s duty of care

- The operating company is responsible for ensuring the suitability of the disinfection method for the item being disinfected (compatibility and intended purpose).
- In particular, the operating company must ensure that the device is only ever used as intended. If the device is used in a way that deviates from its intended use, this can lead to damage or risk for which the manufacturer assumes no liability.
- In order to uphold the function and safety warranty, only original parts from the manufacturer may be used where required.
- The device must only be operated by sufficiently qualified and authorised personnel. These persons should be trained in all applicable aspects of occupational safety and environmental protection on a regular basis. They must be familiar with the instructions for use, and in particular with the safety instructions contained within it.
- The PlasmaEgg® must only be operated when in perfect, fully functional condition.

4. Structure and components

4.1 Control panel

- a) Stand-by button (for switching the device on and off)
- b) Status LED (displays the programme progress – time counts down)
- c) LOG button (for selecting different programme running times)
- d) START button (for starting the selected programme)

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5. Transport and installation

5.1 Transport

The PlasmaEgg® has an integrated carrying handle on the back

5.2 Installation

- Install on a flat surface. Device must be positioned so that the plug can be easily reached.
- Insert the supplied plasma cartridge before first use (see chapter 7).
- Plug in the device.
- Install the device in a side room with a room volume of at least 36 m³.
- Using the PlasmaEgg® for disinfection generates ozone, which is detectable in extremely low concentrations. The device has an integrated filter system that neutralises the plasma.
- Tests to verify compliance with recommended limit values for occupational safety have been carried out and passed.

6. Disinfection with the PlasmaEgg®

6.1 Preparation

- Place the item for disinfection into the chamber. Please make sure that any heavy soiling is removed from items before these are disinfected in the PlasmaEgg®.
- The items must be dry.
- Items can be placed in the chamber in any position/alignment.
- All surfaces that are exposed to ambient air can be disinfected.
- Porous materials are permeated by cold plasma.

6.2 Switching on the device

To switch on the device, press the stand-by button . The status LED changes from green to white. The device is now ready to use. The most recently selected programme is displayed.

6.3 Selecting a programme

The LOG button can be used to select the desired disinfection programme.

There are 3 programmes to choose from (25, 65 and 95 minutes).

Press the LOG button again to switch to the next programme. The suitable programme running time can be found in Section 2.1 (“Performance characteristics”).

6.4 Starting the programme

- Press down the lid, making sure that it is securely closed.
- The disinfection process is started by pressing the START button.
- The status LED flashes purple while the disinfection process is being performed.
- The progress of the programme can be seen at all times via the time display.
At the end of the disinfection process, the display shows a remaining time of 0 minutes, and the status LED changes from purple to white.

6.5 Removing the disinfected item


As soon as the status LED turns white, the lid can be opened and the item removed. Make sure to observe appropriate hand hygiene when removing the disinfected item.

6.6 Switching off the device

The PlasmaEgg® is switched off by pressing the stand-by button . The status LED lights up green. We recommend disconnecting the device from the mains if it is not going to be used for an extended period.

7. Replacing consumables

- The device has been designed such that minimal cleaning and care is required. It is maintenance-free.
- The plasma cartridge needs to be replaced after approximately 250 operating hours. This is indicated by the code E03 on the display. The device does not start anymore.
- We recommend keeping a spare plasma cartridge in stock.
- The plasma cartridge can be replaced by the user. Disconnect device from the mains. Reach under the plasma cartridge and pull it upwards. Push the new cartridge down as far as it will go until it is flush.

8. Hygiene

- Always disconnect the mains plug before cleaning and/or disinfecting the device.

 **Do not spray or pour any liquid directly onto or into the device.**

8.1 Cleaning

- The housing and chamber can be cleaned using a damp cloth and soap or a mild detergent.
- After cleaning, rub with a soft, dry cloth. The device must be dry before it can be used again.

 **Never rinse or flush out the device!**

8.2 Disinfection

- Use alcohol-based disinfectant wipes to disinfect the housing and the chamber.
- Then rub with a clean cloth. The device must be dry before it can be used again.

9. Conditions for storage and transport and ambient conditions

- Prevent all types of shock and impact, regardless of whether or not the device is packaged up.
- Prevent drastic temperature changes, high humidity levels or dusty environments.
- The device can be operated at an ambient temperature of between +5 °C and +40 °C, a relative humidity of between 15% and 90% (non-condensing) and an air pressure of between 700 and 1060 hPa. Do not use the device above 2000 meters above sea level.
- Store and transport the device at an ambient temperature of -40 °C to +70 °C, a maximum relative humidity of 85% condensing and an air pressure of between 500 hPa and 1060 hPa.
- Only operate the device in a dry environment.

10. Warranty

- A 24-month warranty is provided from the date of purchase.
- The warranty is rendered invalid by improper use, if the device is opened, or if non-approved accessories or spare parts are used.
- Signs of use such as scratches and discolourations do not constitute reasons for a warranty claim.

11. Liability

No liability is accepted for damage to the device and other items caused by operating errors or due to a failure to observe the instructions for use. Modifications to the device – in particular internal technical modifications by non-authorised persons – will completely void any warranty claims and will invalidate the product liability.

12. Troubleshooting (faults and remedies)

Code/fault	Possible cause	Remedy
E01	Lid incorrectly attached	Re-position the lid. If the sealing ring is stiff, apply talcum powder.
E02	Plasma source incorrectly inserted	Re-insert until the cartridge sits flush
E03	Plasma cartridge has reached the maximum number of operating hours	Replace plasma cartridge
	Plasma cartridge defective	Replace plasma cartridge
E04	Fan not working	Contact customer services
Device not starting up (status LED not illuminated)	Mains cable not plugged in	Check the plug connection between the mains power and the device. Status LED should light up green when connected to the mains.
	Internal power supply defective	Contact customer services

- If these suggestions do not solve your problem, contact your dealer or the manufacturer.
- The device is maintenance-free. It will need to be sent to the manufacturer for any repairs.
- The regular electrical safety check (recommended every 2 years) can be carried out either by the manufacturer or another authorised person with the necessary specialist knowledge and the appropriate measuring and testing equipment.

The device may only be opened by the manufacturer.

13. Environmental protection and disposal

Before it can be disposed of, the device must be decontaminated and cleaned for the protection of people, the environment and materials. At the end of its product service life, please dispose of the device according to your local waste regulations for electrical appliances.

The device must not be disposed of as household waste.

14. Technical specifications

The expected operating service life is 5 years.

PlasmaEgg®

Rated voltage:	100-240 V AC -15/+10% 50/60 Hz max. 0.6-0.3 Arms
Power consumption:	max. 18 W
Dimensions:	19.5 x 19.5 x 24.5 cm (L x W x H)
Chamber dimensions:	10.0 x 15.3 x 16.0 cm (L x W x H)
Chamber volume:	1.64 l / 1640 cm ³
Maximum fill weight:	10 kg
Programmes:	3 programmes (25, 65, 95 minutes)

Weight:	1.5 Kg
Mains cable:	2-Wire 2 x 0.75 mm; H03VVH2-F; Length: 1.8 m with Euro plug (REF 50083007)
Overtoltage category:	II
Protection class (IP code):	IP 21
Pollution degree:	PD 2
Short-term exposure level:	0,101 ppm
Long-term exposure level:	0,000 ppm

15. Electromagnetic compatibility (EMC)

The PlasmaEgg® device described in the instructions for use must only be operated in accordance with the following information and using the accessories and spare parts specified in Section 17 and only in an electromagnetic environment that corresponds to a domestic environment (house, doctor's practice...) and that has access to the public supply grid. The PlasmaEgg® can also be used in the environment of a professional healthcare facility (medical centre).

It is not permitted to operate the PlasmaEgg® in modes of transport on land, on water and in the air, in specific HF screening areas (e.g. MRT) where there are high-intensity electromagnetic interference levels, or in the direct vicinity of active HF surgical devices or short-wave therapy devices.

WARNING: The use of this device immediately next to other devices or stacked with other devices should be avoided, as this could lead to defective operation. If this type of use is necessary, however, this device and other devices should be observed in order to ensure that they are working correctly.

WARNING: the use of accessories, converters and cables other than those specified or provided for this device by Kimetec may lead to increased electromagnetic emitted interference or reduced electromagnetic interference immunity on the part of the device.

WARNING: Mobile HF communication devices (radio devices) including accessories such as antenna cables and external antennas should not be used at less than 30 cm from the PlasmaEgg® parts and cables designated by Kimetec.

Failure to observe this rule may reduce the performance features of the device.

15.1 Behaviour of the PlasmaEgg® in the described electromagnetic environment

If the device is operated in the electromagnetic environment described below and with the accessories/spare parts specified in the instructions for use, there are no restrictions in terms of its key performance features and no impairment of its basic safety. Nevertheless, electrostatic discharge may in certain circumstances mean that a slow or non-responsive display may be expected. The device can be restarted by pressing the stand-by button. Leave at least 5 seconds between switching off and switching back on again.

15.2 Description of the electromagnetic environment

The device is intended for operation in the electromagnetic environment described below. The operating company or user is responsible for ensuring that the device is only operated in this environment.

15.3 Emissions

Phenomenon/measurements	Conformity	Electromagnetic environment/guideline
Conducted and radiated emissions CISPR 11	Group 1	The device is designed for use in all facilities, including residential homes, that are connected directly to the public mains that also supplies buildings used for residential purposes, as well as for use in professional healthcare facilities
	Class B	
Harmonic distortion IEC 61000-3-2	Class A	
Voltage fluctuations and flicker IEC 61000-3-3	Passed	

The device uses HF energy exclusively for its internal function. Its HF emissions are therefore very low, and it is unlikely that nearby electronic equipment will suffer any interference, with the exception that the frequency range of 13.56 Mhz is modulated with a bandwidth of 1.7 Mhz ASK. The transmitter has a radiation power of 11.5 mW corresponding to 10.6 dBm (Rout 20 Ohm;RL 50 Ohm U0pp = 3 V).

15.4 Immunity

EMC basic standard/ phenomenon	Interference immunity test level	Electromagnetic environment/guideline
IEC 61000-4-4 Electrical fast transients (burst)	± 2 kV contact 100 kHz repetition frequency	The device is designed for use in facilities that are connected directly to the public mains that also supplies buildings used for residential purposes, as well as for use in professional healthcare facilities
IEC 61000-4-5 Surges (Surge)	± 0.5 kV, ± 1 kV LgL LgE omitted due to SKII	The device is designed for use in facilities that are connected directly to the public mains that also supplies buildings used for residential purposes, as well as for use in professional healthcare facilities
IEC 61000-4-6 Conducted disturbances, induced by high-frequency fields	3 V 0.15 MHz to 80 MHz 6 V in ISM bands 6.765 MHz - 6.795 MHz 13.553 MHz - 13.567 Mhz 26.957 MHz - 27.283 MHz 40.66 MHz - 40.70 MHz and 6 V in amateur radio bands 1.8 MHz - 2.0 MHz 3.5 MHz - 4.0 Mhz 5.3 MHz - 5.4 MHz 7 MHz - 7.3 MHz 10.1 MHz - 10.15 MHz 14 MHz - 14.2 MHz 18.07 MHz - 18.17 MHz 21.0 MHz - 21.4 Mhz 24.89 MHz - 24.99 MHz 28.0 MHz - 29.7 MHz 50.0 MHz - 54.03 MHz 80 % AM at 1 kHz	The environment can correspond to a domestic environment or a professional healthcare facility.
IEC 61000-4-11 Voltage interruptions (smallest and largest rated voltage)	0 % UT; 250 / 300 periods	If the user requires continued operation even in the event of interruptions to the power supply, it is recommended that the device be powered via an interruption-free power supply (UPS) or battery.

EMC basic standard/phenomenon	Interference immunity test level	Electromagnetic environment/guideline
IEC 61000-4-1 Voltage drops (smallest and largest rated voltage)	0 % UT; 1 / 2 periods At 0, 45, 90, 135, 180, 225, 270 and 315 degrees	The device is designed for use in facilities that are connected directly to the public mains that also supplies buildings used for residential purposes, as well as for use in professional healthcare facilities.
	0 % UT; 1 period and 70 % UT; 25 / 30	
	0 % UT; 1 period and 70 % UT; 25 / 30 periods Single-phase: at 0 degrees	
IEC 61000-4-3 High-frequency electromagnetic fields in the direct vicinity of wireless communication equipment	10 V/m 80 Mhz to 2.7 Ghz 80 % 1 khz	The environment can correspond to a domestic environment or a professional healthcare facility.

Safety distances to certain wireless communication devices

EMC basic standard		Phenomenon				
IEC 61000-4-3		High-frequency electromagnetic fields in the direct vicinity of wireless communication devices				
Test frequency MHz	Frequency band MHz	Radio service	Modulation	Maximum output W	Distance m	Interference immunity test level V/m
385	380 - 390	TETRA 400	Pulse modulation with square wave signal dt 50 % 18 Hz	1.8	0.3	27
450	430 - 470	GMRS 460; FRS 460	Frequency modulated ± 5 kHz swing 1 kHz sine	2	0.3	28

EMC basic standard		Phenomenon				
IEC 61000-4-3		High-frequency electromagnetic fields in the direct vicinity of wireless communication devices				
Test frequency MHz	Frequency band MHz	Radio service	Modulation	Maximum output W	Distance m	Interference immunity test level V/m
710	704-787	LTE band 13, 17	Pulse modulation with square wave signal dt 50 % 217 Hz	0.2	0.3	9
745						
780						
810	170-1990	GSM 1800, GSM 1900 CDMA 1900, DECT LTE band 1, 3, 4, 25; UMTS	Pulse modulation with square wave signal dt 50 % 217 Hz	2	0.3	28
870						
930						

2450	2400–2570	Bluetooth, WLAN 802.11 b/g/n, RFID 2450 LTE band 7	Pulse modulation with square wave signal dt 50 % 217 Hz	2	0.3	28
5240	5100–5800	WLAN 802.11 a/n	Pulse modulation with square wave signal dt 50 % 217 Hz	0.2	0.3	9
5500						
5785						

EMC basic standard/ phenomenon	Interference immunity test level	Electromagnetic environment/guideline
IEC 61000-4-8 Power frequency magnetic fields	30 A/m	No impairment
IEC 61000-4-39 Magnetic fields in close range	30 kHz CW 8A/m 134.2kHz PWM 2.1kHz 65A/m 13.56 Mhz PWM 50Khz 7.5A/m	No impairment
IEC 61000-4-2 Static electricity discharge	8 kV contact ± 2 kV, ± 4 kV, 8 kV, ± 15 kV air	Floors must be made of wood or concrete, or be fitted with ceramic tiling. If the flooring is fitted with synthetic material, the relative humidity must be at least 30%.

16. Legend

Symbols on the packaging, the type label and in the instructions for use

	Safety information		Warning information / Caution
	Prohibited action symbol		Manufacturer
	Store in dry conditions		Protect against sunlight
	Do not dispose of as household waste		Applicationpart type BF
CE	Conformity symbol		Do not open
	Device in protection class II	IP21	Protected against the ingress of solid foreign bodies ≥12.5 mm diameter. Protected against vertically falling drops.
	Item number		Serial number
	Follow the instructions for use		Date of manufacture

Symbols on the control panel

	Stand-by button		START and LOG button
---	-----------------	---	----------------------

17. Accessories and spare parts

Only use accessories and spare parts that meet the manufacturer's specifications.

Designation	REF
PlasmaEgg®	500610
Mains cable	50083007
Lid	500611
Plasma cartridge	500612

Kimetec GmbH
Gerlinger Str. 36-38
71254 Ditzingen
Germany

Tel.: +49 (0) 71 56 / 1 76 02 - 200
Fax: +49 (0) 71 56 / 1 76 02 - 500
info@kimetec.de