

The colorful single-use tourniquet for children

Optimal compression force (approx. 20 mmHg)

- By extension of the elastic strap, the optimal compression force is reliably achieved
- Defined compression force – supporting the user
- Complications of venipuncture (e.g. rupture) can be reduced

No laborious cleaning or disinfection

- Single-use product - no cleaning or disinfection necessary
- No cross-contamination with MRSA and other germs

Comfort for medical professionals

- Easy handling even with gloves
- Micro velcro for easy closing and opening
- Particularly low weight (2.5 g/pc.)

Comfort for the patient


- Velvety soft fleece – no bruising, lacing or constricting of the skin
- Free from latex, silicone and plasticisers

Customized and colorful

- Starting from 20,000 pieces – customized colors and designs


Areas of use

- Isolation wards
- House calls
- Paediatrics
- Oncology




The colorful single-use tourniquet for children

Application



Place **unprinted side** on the patient's arm.



Build up compression force by extending the tourniquet completely.



Maintain the tension!
Place tightly around the arm and attach in the velcroed area.



Simply pull to release the strap.

Do not reuse to avoid cross-contamination and dispose of in accordance with regulations.

Video



Ordering data



Medical device

Product description:	CBC VECOS® KIDS
Strap	3 cm wide, 50 cm long
REF 521003 013	125 pieces in dispenser box
REF 521003 015	1 poly bag (refill pack) with 125 pieces
REF 521003 014	8 boxes with 125 pieces in each dispenser box
REF 521003 016	8 x 1 poly bag (refill pack) with 125 pieces

Save costs with the CBC VECOS® KIDS:

- No reprocessing, cleaning, or disinfection.
- Free from latex, silicone, or DEHP (plasticizer).