

## Desinfektionsmittelliste für CBC NUCOS

### List of disinfection liquids for CBC NUCOS

Verträgliche Desinfektionstücher: <i>compatible disinfection wipes:</i>	Verträgliche Flächendesinfektionsmittel: <i>compatible disinfection fluids:</i>
<ul style="list-style-type: none"><li>• Ecolab Incidin™ Alcohol Wipe</li><li>• Hartmann Bacillof® Sensitive Tissues</li><li>• Schülke mikrozyd® universal wipes premium</li><li>• Bode Mikrobac® Tissues</li><li>• Ecolab Incidin™ OxyWipe S</li><li>• Dr. Schumacher Cleanisept® Wipes Maxi</li><li>• Dr. Schumacher Ultrasol OXY® Wipes</li></ul>	<ul style="list-style-type: none"><li>• Bode Kohrsolin® FF</li><li>• Ecolab Incidin™ Rapid</li><li>• Ecolab Incidin™ Active (powder)</li><li>• Ecolab Incidin™ Pro</li><li>• Schülke perform® (powder)</li></ul>

Unverträgliches Flächendesinfektionsmittel: <i>Non-compatible disinfection fluid:</i>
<ul style="list-style-type: none"><li>• Bode Bacillof® AF</li></ul>