


CBC Slingcross®

Single-use tourniquet for children The basic product for venipuncture


Advantages at a glance:

- Improves patient hygiene
- Single use of strap:
no cross-contamination with MRSA and other germs
- Allows individual pressure for every patient
- The velvety soft fleece is gentle on the skin
- Easy handling even with gloves
- Micro velcro for easy closing and opening,
sticks along the complete strap
- Also usable as compress after venipuncture
- Convenient strap removal from the roll
- **Made in Germany**
- Free from latex, silicone and plasticiser



Areas of use:


- Pediatrics
- Isolation wards
- Rehab centres
- Home calls
- Oncology, etc.




CBC Slingcross®

Single-use tourniquet for children The basic product for venipuncture


Application:



Place **with the velcro outside** on the patient's arm.



Tighten the strap – pulse needs to be traceable – and attach in the velcroed area.



Simply pull to release the strap.

Dispose of single-use tourniquet after its use in accordance with regulations to avoid cross-contamination.



Ordering data:



Medical device

Product description	CBC Slingcross® children
Strap	1,5 cm wide, 30 cm long
REF 521003 007	1 roll with 50 pcs. in dispenser box

Save costs with CBC Slingcross®:

- No reprocessing, cleaning, or disinfection.
- Additional benefit: serves as plaster and bandage.
- Allergy to latex, silicone, or DEHP excluded.